

## “OffBead” Protocol for Capturing in Solution

**NOTE: Before proceeding this protocol, please read carefully the instruction in the handbook “Guideline for CCMS”! This is only a brief description of the procedure!**

**Perform steps 1 - 8 at 0 - 4 °C using the caproBox™ for cooling.** Avoid bubbles during handling which may cause denaturation of proteins. We recommend 2 - 5 mg/ml total protein concentration for cell lysates. The assay is designed to use either single 200 µl PCR tubes or 200 µl PCR tube strips.

- 1) Mix the sample with all reaction components in a total volume of 100 µl as described below.



**NOTE:** The schemes below are examples to prepare the reactions. For your own experiments, please use the provided pipetting scheme (see below) or the reaction volume calculator (see at: [www.caprotec.com/support/downloads](http://www.caprotec.com/support/downloads)).



**NOTE:** For the efficiency of the CCMS approach, place the PCR tube strips in the caproBox for cooling, add all components **exactly** following the sequence below - mix gently after each addition.

### **Example to prepare Capture Compound assay (“A”) and Competition control (“C”)**

Component	Stock conc.	Assay conc.	Capture Compound assay “ <u>A</u> ”	Competition control “ <u>C</u> ”
H <sub>2</sub> O			39.3 µl	14.3 µl
cAMP competitor	4 mM	1 mM	--	25 µl
5x CB	5 x	1 x	20 µl	20 µl
PKA RI*	16 µM	0.11 µM	0.7 µl	0.7 µl
Cell lysate	10 mg/ml	3 mg/ml	30 µl	30 µl
<b>10 min pre-incubation</b> <b>Proceed to add:</b>				
cAMP-CC	100 µM	10 µM	10 µl	10 µl
<b>Total volume</b>			100 µl	100 µl

\*For the initial experiments, it is recommended to use the **positive control** (PKA RI) in addition to the Capture Compound™ assay “A” and the Control “C”. **Additional control reactions** are listed under “*checklist for OffBead protocol*”. In addition, keep a 1 µl sample of the prepared Capture Compound assay “A” solution as a reference for subsequent SDS-PAGE and MS analysis.

- 2) Incubate solutions “A”, “C” and in case additional control reactions (cf. “Checklist”) on a rotating wheel at 4 °C for 30 min (or gently re-suspend every 5 min).

- 3) Irradiate for 10 min using the caproBox. *Caution: Lids of PCR tube strips could be warm! If necessary, cool down the lids of PCR tube strips before proceed step 4 to avoid protein precipitation.*



**Note:** Before irradiation, make sure that no suspension is in the lids of PCR tubes.

- 4) Add 25 µl cAMP competitor solution to “**A**”, gently mix and incubate for 10 min.



**Note:** Adding free competitor after photo cross-linking will displace all non-covalently bound proteins. Covalently bound proteins will not be affected by free competitor. This step can be omitted if the Capture Compound experiment is paired with a pull down assay: All affinity driven, non-covalently bound proteins will be isolated in addition to the covalently attached proteins after UV-activation of the reactivity function (cf. “Checklist”).

- 5) Add 25 µl 5x WB1 to “**A**”, “**C**” and additional control reactions and gently homogenize.
- 6) Add 50 µl well resuspended Streptavidin magnetic beads (SA-MB) to “**A**”, “**C**” and any additional control reactions. Gently mix the resulting suspensions.
- 7) Incubate for 30 min on a rotating wheel at 4 °C (or gently re-suspend every 5 min).



**Note:** SA-MB must stay in suspension.


- 8) Collect SA-MB protein complexes (beads) from the reactions by using caproMag™. Allow collection process to proceed for 2 min, until supernatant appears clear. Discard tubes containing supernatant. Add 200 µl WB1 in new PCR tubes and gently resuspend the beads.



**Note:** Take care not to pinch your fingers between the Neodymium magnet and the steel plate. Keep away from pace makers or other metallic objects. For further handling advices, please download “cartoon for using the caproMag” from [www.caprotec.com/support/downloads](http://www.caprotec.com/support/downloads)

- 9) Collect the beads with the caproMag and repeat the washing step five times without changing the reaction tube by vigorously mixing the reaction suspensions as defined in step 8. Discard supernatants after all wash steps. Wash the beads once with 200 µl ultrapure water.

The collected beads in step 9 are now ready for further analysis. For sample storage add 10 µl ultrapure water and keep at 4 °C for up to one week.



### Checklist for “OffBead” Protocol and Recommended Additional Control Reactions

	<u>A</u>	<u>C</u>	PD	C-PD	<u>A</u> + PD	C1-MB	C2-MB	C3-MB	C4-MB	C5-MB
H <sub>2</sub> O	✓	✓	✓	✓	✓	-	✓	✓	✓	✓
Competitor (prior irr.)	-	✓	-	✓	-	-	-	-	-	-
5x CB	✓	✓	✓	✓	✓	-	✓	✓	✓	✓
PKA RI*	✓	✓	✓	✓	✓	-	✓	✓	✓	✓
Cell lysate	✓	✓	✓	✓	✓	-	✓	-	✓	-
cAMP-CC	✓	✓	✓	✓	✓	-	-	-	-	-
Biotin	-	-	-	-	-	-	-	✓	-	✓
Irradiation	✓	✓	-	-	✓	-	-	-	✓	✓
Competitor (after irr.)	✓	-	-	-	-	-	-	-	-	-
5x WB1	✓	✓	✓	✓	✓	-	✓	✓	✓	✓
SA-MB	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Incubation	✓	✓	✓	✓	✓	-	✓	✓	✓	✓
Wash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

CB = capture buffer; cAMP-CC = cAMP Capture Compound; WB1 = wash buffer 1; SA-MB = Streptavidin magnetic beads;

\*For the initial experiments, it is recommended to use a positive control (PKA RI) in addition to the Capture Compound assay “A” and the Control “C”.

Synonym	Description	Comments
<u>A</u>	Capture Compound assay	
<u>C</u>	Control of Capture Compound assay	Proteins only detected in <u>A</u> or detected in a much lower extent in <u>C</u> are target binders. The proteins are covalently cross-linked to CC. The cross-link position within the protein sequence may be determined by MS.
PD	Pull down assay	
C-PD	Control of pull down assay	Proteins only detected in PD in comparison to C-PD are strong or highly abundant target binders. They are not removed from the SA-MB by washing steps. No covalent cross-link between CC and proteins occurred.
<u>A</u> + PD	Combined Capture Compound and pull down assay	Use <u>C</u> as control.
C1-MB	Control, examining Streptavidin magnetic beads (SA-MB) only.	Only Streptavidin should be detected.
C2-MB	Control, examining CC-independent binding to the SA-MB	Biotinylated proteins are detected if not removed by pre-treatment of the lysate. Additionally, proteins with unspecific binding to the SA-MB may be detected (compare C3-MB).
C3-MB		Control, examining CC-independent binding to the SA-MB excluding biotinylated proteins
C4-MB		Control, examining CC-independent binding to the SA-MB including irradiation induced unspecific binding
C5-MB		Control, examining CC-independent binding to the SA-MB including irradiation induced unspecific binding and excluding biotinylated proteins

## Contact

### Headquarters

**caprotec bioanalytics GmbH**  
Volmerstrasse 5  
D-12489 Berlin

Phone: +49 30 63 92 39 90

Fax: +49 30 63 92 39 89

Web: [www.caprotec.com](http://www.caprotec.com)

Email: [sales@caprotec.com](mailto:sales@caprotec.com)

**caprotec Inc., USA**

15 New England Executive Office Park  
Burlington, MA 01803, USA

Phone: + 1 781 685 4992

Fax: +1 781 685 4601

Web: [www.caprotec.com](http://www.caprotec.com)

E-Mail: [info@caprotec.com](mailto:info@caprotec.com)

### Copyright

© 2008 caprotec bioanalytics GmbH. All rights reserved. Reproduction in whole or in part only with permission of caprotec bioanalytics GmbH.

### Trademarks

caprotec, caproKit, Capture Compound, caproBeads, caproBox, caproMag and ImproMed are registered trademarks of caprotec bioanalytics GmbH. All other used tradenames or trademarks belong to their respective proprietaries.

**FOR RESEARCH ONLY. NOT FOR HUMAN OR DIAGNOSTIC USE.**

### Products & Services

CCMS technology is made available as ready to use caproKit reagents and services.

For more information please visit [www.caprotec.com](http://www.caprotec.com)



Or contact us. Email: [info@caprotec.com](mailto:info@caprotec.com)

Phone: +49 30 6392 4004

Products and Services are for Research use only.

© 2008 caprotec bioanalytics GmbH

# “OffBead” Protocol for Capturing in Solution



## Pipetting Scheme - OffBead

### Preparing Capture Compound assay ("A") and Competition control ("C")

Component	Stock conc.	Assay conc.	Capture Compound assay " <u>A</u> "	Competition Control " <u>C</u> "
H <sub>2</sub> O				
cAMP competitor	4 mM	1 mM	--	25 µl
5x CB	5 x	1 x	20 µl	20 µl
PKA RI*	16 µM	0.11 µM	0.7 µl	0.7 µl
Cell lysate				
<b>10 min pre-incubation</b> <b>Proceed to add:</b>				
cAMP-CC	100 µM	10 µM	10 µl	10 µl
<b>Total volume</b>			100 µl	100 µl

\*For the initial experiments, it is recommended to use the **positive control** (e.g. PKA RI) in addition to the Capture Compound assay "A" and the Control "C".